

INSTRUCTIONS FOR USE

Please read the instructions carefully before performing the test. Follow the instructions, and do not modify the process. Strict adherence to the guidelines will avoid inaccurate results and achieve optimal performance of Saligen.

Product name

Saligen

Intended use


Saligen is an in vitro diagnostic medical device based on the Immuno-chromatographic assay (ICA) principle for the qualitative detection of SARS-CoV-2 antigens in human saliva. This test is used to detect antigens of the SARS-CoV-2 virus in people suspected of COVID-19. This product is intended exclusively for professional use in the laboratory or at the point-of-care.

Summary and explanation of the test

COVID-19 is a respiratory disease caused by a new type of coronavirus (SARS-CoV-2) first identified in December 2019 in Wuhan, China. Common signs of infection include, but are not limited to, respiratory symptoms, fever, cough, shortness of breath, reduced sense of smell or taste. In severe cases, the infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and death. Coronaviruses are a group of viruses that cause symptoms from the common cold to more severe illnesses such as Middle East Respiratory Syndrome (caused by MERS-CoV) and Severe Acute Respiratory Syndrome (caused by SARS-CoV).

Principle of the procedure

Saligen uses COVID19 antibodies, which are labeled with small gold particles and are attached to a nitrocellulose membrane near the sample hole of the test card (see also illustration below). After its application, capillary forces are pulling the sample from the sample hole to the test region of the device. When the liquid of the sample reaches the COVID19 antibodies, they detach from the membrane and are moved along the test card.



If the sample contains SARS-CoV-2 antigens ("analyte"), these bind to the labelled antibodies to form analyte-labeled antibody complexes. When these complexes reach the test line of the test card, they are retained on the test line by another set of COVID19 antibodies, which are immobilized on the nitrocellulose membrane. These so-called sandwich complexes appear as a color band on the test line. If the sample does not contain SARS-CoV-2 antigens, no sandwich complexes are formed and no color band appears on the test line. Regardless of the presence or absence of SARS-CoV-2 antigens in the sample, a color band will appear on the control line of the test card. If no color band appears on the control line, the test card has not worked as intended.

Kit Components

- 1) Test card
- 2) Extraction buffer tube
- 3) Filter cap
- 4) Specimen cap collector

Required materials not included

- Timer or stopwatch

Kit storage and stability


- Saligen should be stored at 2-30°C in a dry place. When stored and handled as directed, the test cards an reagents are stable until expiration date indicated on kit labels.
- Test cards should be used immediately after opening the pouch.

Sample collection

Before the collection, do not eat, smoke, chew or drink any beverages apart from water.

Saliva specimens

- The person to be tested collects saliva in the mouth on the tip of the tongue for 30 seconds (approximately 0.5 mL); see also illustration below.
- Spit the collected saliva into the extraction buffer tube directly for immediate use. The applicator can be used to assist this step. By adding the saliva, the volume in the tube should approximately double.



- Do not use stored specimens. Long-term storage may result in a signal decrease.
- Do not freeze the sample. Multiple freeze/ thaw cycles may result in a signal decrease.

Warnings and precautions


- This product is intended for in vitro diagnostic use.
- This product is intended for single use.
- This product is intended for professional use.
- This product is intended for POCT use with human saliva
- Assay should be performed as directed in the instructions for use to obtain accurate results.
- Do not use beyond the expiration date or damaged products.
- Do not use any other reagents that are not provided in this kit and do not mix components of different lots.
- This reagent can be stored at room temperature (15-25°C). Reagents stored or samples collected at lower temperatures should be allowed to come to room temperature before use.
- Remove the test card from the pouch and use it as soon as possible to avoid prolonged exposure to air. Prolonged exposure to air affects the test results.
- Follow laboratory test procedures for infectious diseases. Waste after use should be treated as infectious material and not disposed of randomly.
- Appropriate safety assurance procedures should be in place for infectious agents and materials.
- Wear gloves to handle samples and reagents.
- Do not suck the samples and reagents.
- Do not smoke, eat, drink, use cosmetic or touch contact lenses while handling the product.
- Spilled samples or reagents should be cleaned with disinfectants.
- Disinfect and dispose of all samples, reagents, and potential contaminants following applicable local regulations.

Preparation for use

Reagents should be allowed to stand at room temperature for 20-30 minutes before testing. Do not use samples, which have been stored for prolonged times after collection.


Assay procedure for saliva specimens

1. Collect the sample as directed in the "Sample collection" and "Saliva specimens" section.
2. Cover the tube with a filter cap and tighten the lid. Mix the contents by turning the tube upside down 10 times.
3. Open the test card pouch and place the test card on a flat surface. Apply a few drops of the saliva extraction buffer mix into the sample hole of the test card. The sample hole should be almost completely filled. Make sure not to use less than 2-3 drops.
4. Read the results after 10 - 15 minutes.



Reading the test card later than 20 minutes after applying the sample diluent may give inaccurate results.

Interpretation of results



Using the test card can lead to three different results:

1. If only one color band appears in the test region near the letter "C", the result is valid and "non-reactive", meaning no SARS-CoV-2 antigens could be detected.
2. If a second color band appears in the test region near the letter "T", the result is valid and "reactive", meaning SARS-CoV-2 antigens were detected.
3. If no color band appears or if only one color band appears near the letter "T", the result is invalid. In this case the result cannot be used, because the test did not work as intended. See section "Internal Control" for details.

Internal Control

Saligen test contains a built-in internal control in the test card. A color band appearing in the control region (C) is designed as an internal control. The appearance of the control line confirms that sufficient flow has occurred, and that the test card is working normally. If the control line does not appear within 10 minutes, it is considered an error in the test result and it is recommended to test again with the same sample and a new device. If there is again no color band on the internal control line on the retest, contact the manufacturer or distributor.

External Controls

- External Positive and Negative controls may be used with the test kit. These controls provide additional quality control material to assess that the test reagents react as expected. Positive controls shall lead to "reactive" results and Negative controls shall lead to "non-reactive" results.
- Controls are recommended to be run once for each new kit lot.
- Prepare solutions for Positive controls following the instructions provided with the control material.
- For external Negative controls it is recommended to use the Extraction buffer included in the kit.
- Perform controls using the same procedure as used for patient specimens.
- If the kit controls do not perform as expected, do not report patient results. Contact the manufacturer or distributor.

Limitations of the procedure

- The results Saligen should not be considered as absolute, and shall not be the sole basis for treatment or patient management. The infection should be confirmed by a specialist along with other experimental results, clinical symptoms, epidemiology, and additional clinical data.
- This kit detects both SARS-CoV and SARS-CoV-2, regardless of their viability. This kit does not differentiate between SARS-CoV and SARS-CoV-2.
- In the early stages of infection, low levels of antigen expression can result in non-reactive results.
- Due to the limitation of the assay methods, non-reactive results cannot entirely rule out the possibility of infection.
- This product can qualitatively detect SARS-CoV or SARS-CoV-2 antigens in human saliva and can't not determine the specific antigen quantity in the sample.

Performance characteristics

Limit of detection (LoD)

The LoD was determined using limiting dilutions of inactivated SARS-CoV-2 (ZeproMetrix, #0810587CFHI) in two separate methods. The inactivated virus was spiked into the extraction buffer processed with a non-reactive saliva sample to have a concentration of TCID50 of 1.15 x 10⁶/ml. Each sample was serially 10-fold diluted and by testing in triplicate, a tentative LoD showing 100% (3/3) reactive rate was determined for each. For confirmation LoD study, 4 concentrations below the lowest concentration of the pre-test were tested in 20 replicates and a concentration showing over 100% (20/20) reactive results was determined as the LoD of the Saligen for each.

- Saliva LoD: 1.44 x 10³ TCID50/ml Cross-reactivity/ Microbial interference Viruses/bacteria listed below were confirmed not to have cross reactivity or cause interference with Saligen.
- Virus (10⁵ TCID50/mL): Adenovirus type 1, Adenovirus type 7, Coronavirus 229E, Coronavirus NL63, Coronavirus OC43, MERS-CoV, Cytomegalovirus, Influenza A H3N2, Influenza A H1N1, Influenza B, Enterovirus type 71, Parainfluenza type 1, Parainfluenza type 2, Parainfluenza type 3, Parainfluenza type 4A, Measles virus, Human Metapneumovirus, RSV type A, RSV type B, Rhinovirus, Epstein Barr virus, Mumps virus and Coronavirus HKU1
- Bacteria (10⁶ CFU/mL): B. pertussis, E. coli, H. influenzae, M. catarrhalis, C. pneumoniae, L. pneumophila, M. pneumoniae, M. tuberculosis, N. meningitidis, P. aeruginosa, S. epidermidis, S. pneumoniae, S. pyogenes, S. salivarius and S. aureus

Endogenous interference

Potential interfering substances listed below were confirmed not to have a response with Saligen.

- Mucin (4 mg/mL), Human Blood (2%), 4-Acetamidophenol (10 mg/mL), Acetylsalicylic Acid (20 mg/mL), Chlorpheniramine (5 mg/mL), Diphenhydramine (5 mg/mL), Guaiacol glyceryl ether (20 mg/mL), Oxymetazoline (0.05 mg/mL), Phenylephrine (1 mg/mL), Fexofenadine (500 mg/mL), Amantadine (500 mg/mL), Ribavirin (500 mg/mL), Pseudoephedrine HCl (20 mg/mL), Ibuprofen (10 mg/mL) and Tamiflu (48 mg/mL), Naso GEL (5%), Chloraseptic (1.5 mg/mL), Cromolyn (15%), Zicam (5%), Homeopathic preparations (1:10 dilution), Sore Throat Phenol Spray (15%), Tobramycin (4 µg/mL), Mupirocin (10 mg/mL), Fluticasone Propionate (5%).

Lot number	In vitro diagnostics medical device	Consult instructions for use
Sufficient for n tests	Store at 2-30°C	Do not reuse
Expiration Date	Manufacturer	Caution
European authorized representative	Conformity European	
Bakter Medical s.r.o. Chaloupky 171/33, Komin, 624 00 Brno IČO: 09418806 DIČ: CZ09418806 Tel: +420 725 395 600 www.saligen.com	Bakter Medical s.r.o. Chaloupky 171/33, Komin, 624 00 Brno IČO: 09418806 DIČ: CZ09418806 Tel: +420 725 395 600 www.saligen.com	

Saligen